


Detection & Speciation of Organic Peroxy Radicals (ROO•)



Gas Kinetics 2024
14 July - 18 July 2024, Leeds, UK

Olivier Durif¹, Barbara Nozière¹

¹ Kungliga Tekniska Högskolan, Department of Chemistry, Stockholm, Sweden



Keywords : *Organic Peroxy Radical, Proton-Transfer-Reaction Mass Spectrometry, Flow Tube, Gas-Phase, Chemical Kinetics*

Take home message

A new experimental setup enables the direct and specific measurement of the simplest organic peroxy radicals in the laboratory ($\text{CH}_3\text{OO}^\bullet$, $\text{C}_2\text{H}_5\text{OO}^\bullet$, $\text{C}_3\text{H}_7\text{OO}^\bullet$, ...). Measurements demonstrated a strong uptake of organic peroxy radicals on metal, indicating redox reactions.[1]

Avoid all metal surfaces and use instead a fluoropolymer to maximize ROO• detection!

Research Topic

We have set up an experiment to study kinetically gas-phase reactions involving radicals at room temperature and atmospheric pressure.

The kinetics reactor is simple and innovative (Fig. 1). Our method is precise, straightforward, inexpensive, and makes it possible to study surface uptakes. It can also be set up without any movable parts.

In the principle, the flow tube is here an elementary 1/4" tube. The radicals are produced continuously in a UV-C cell by the photolysis of a precursor. The radicals are then flown through a tube before their detection by proton-transfer-reaction mass spectrometry (Fig. 2). It is the residence time in the tube, by changing its length or the flows, that allows for kinetic measurements. Our laboratory approach has recently been applied to measure the uptake of small organic peroxy radicals on metallic and non-metallic surfaces (Fig. 3).

Mass Spectrometry Data




FIGURE 2 – Mass spectra acquired every second allows identification of surface reactions.

Challenges, Present and Near Future

- ▶ Raw data and processing tools available on Zenodo: zenodo.org/records/8250035
- ▶ Complex data processing, GBytes of data acquired every experiment hour. A Julia package has been developed for data processing: github.com/MassSpec.jl
- ▶ Impact of relative humidity on ROO• signal ill-understood. Ongoing studies.
- ▶ An open possibility to investigate ROO• uptake by organic and mineral surfaces.

Experimental Setup




FIGURE 1 – A photolysis cell, a simple tube and a Proton-Transfer-Reaction Mass Spectrometer.

The Uptake Kinetics

Chemical system: $\text{ROO}^\bullet \xrightarrow{\text{Surface}} \text{P}$, $2\text{ROO}^\bullet \xrightarrow{\text{Gas}} \text{P}'$.

When gas-phase matter (Glass and PFA), mixte-order reaction:

$$\frac{d[\text{ROO}^\bullet]}{dt} = -2k_{\text{self}}[\text{ROO}^\bullet]^2 - k'[\text{ROO}^\bullet].$$

Analytical solution: $[\text{ROO}^\bullet](t) = \frac{k' \exp(-k't)}{\frac{k'}{[\text{ROO}^\bullet]_0} + 2k_{\text{self}}(1 - \exp(-k't))}$.

On metals (calibration-free): $[\text{ROO}^\bullet](t) = [\text{ROO}^\bullet]_0 \exp(-k't)$.




FIGURE 3 – Kinetics results demonstrated strong uptakes of organic peroxy radicals on metal, following the electrochemical series.

[1] Olivier DURIF, Felix PIEL, Armin WISTHALER et Barbara NOZIÈRE. « Strong Uptake of Gas-Phase Organic Peroxy Radicals (ROO•) by Solid Surfaces Driven by Redox Reactions ». In : *JACS Au* (avr. 2024). DOI : 10.1021/jacsau.4c00060.

ACKNOWLEDGMENTS. This work is part of the ERC Advanced Grant Project EPHEMERAL (grant agreement No 884532) and has received funding from the European Research Council (ERC) under the European Union's Horizon 2020 research and innovation programme.